



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bode rack (PP) for 18 tubes 26–30 mm, for 2–5166

€91.70
plus VAT &
Shipping

Product Images

Description

Rack made of PP

for 18 tubes of 26–30 mm Ø

Accessories for: Mini incubator (Order No. 2-5165)

Additional Information

No.	2-5166
Manufacturer (Brand)	Bode
EAN	4058072004408
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Material	Polypropylene (PP)
suitable for	Test tubes/tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	30 mm
Capacity	18
Squares (WxL)	KA
Type frame	Test tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes
suitable for cooling in	Freezer/ freezer

