



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic test tube rack for 2 x 6 tubes, 145 mm high

Special Price
€29.00 was
€43.50
plus VAT &
Shipping

Product Images



Description

Glued and screwed, Heidelberg Clinic model, heavy-duty design.

Height: 145 mm

Hole width: 18 mm, two-stage

Additional Information

No.	1-5121
Manufacturer (Brand)	neoLabLine
EAN	4058072035761
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	155 mm
Height	145 mm
Depth	64 mm
suitable for	Test tubes/tubes
Racks features	incubator suitable for different vessel types/volumes
Number of floors	1
DM Bore	18 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 6
Type frame	Test tube rack
for vessels with DM up to	18.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

