



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic test tube rack for 2 x 10 tubes, 145 mm high

**€73.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Glued and screwed, Heidelberg Clinic model, heavy-duty design.

Height: 145 mm

Hole width: 18 mm, two-stage

## Additional Information

---

No.	1-5122
Manufacturer (Brand)	neoLabLine
EAN	4058072035778
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	272 mm
Height	145 mm
Depth	64 mm
suitable for	Test tubes/tubes
Racks features	incubator suitable for different vessel types/volumes
Number of floors	1
DM Bore	18 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 10
Type frame	Test tube rack
for vessels with DM up to	18.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

