




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® test tube rack 5 x 10 places

**€40.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Space for many test tubes!

Test tubes, bottles or plastic tubes can be stored vertically or inverted to be placed on the rods for drying.

Material: Polypropylene

Autoclavable up to 121°C

## Additional Information

---

No.	6-2352
Manufacturer (Brand)	neolab
EAN	4058072025779
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	180 mm
Height	67 mm
Depth	100 mm
suitable for	Test tubes/tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	17 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	5 x 10
Type frame	Test tube rack
for vessels with DM up to	17.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

