



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® test tube rack blue, ABS plastic, Z-
shape, Ø: 17 mm, 225 x 114 x 60 mm, Temp.:
...90°C**

€8.95 As low as
€8.25
**plus VAT &
Shipping**

Product Images



Description

The Z-shaped blue test tube racks are made of robust ABS plastic.

Due to their practical shape, the racks can be plugged together in a particularly space-saving way, when they are not needed.

This leaves you enough space for other applications in your laboratory.

- Ø = 17 mm
- 50 places
- 225 x 114 x 60 mm (L x W x H)
- Temperature resistant up to 90 °C

Additional Information

No.	2-6996
Manufacturer (Brand)	neoLab
EAN	4058072080723
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Acrylonitrile butadiene styrene (ABS)
Wide	225 mm
Height	60 mm
Depth	114 mm
Racks features	stackable/connectable for different vessel types/volumes
Number of floors	1
DM Bore	17 mm
for vessels with volume up to	15 ml
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Squares (WxL)	5 x 10
Type frame	Test tube rack
for vessels with DM up to	17.0 mm
for micro vessel/tube shape	for all shapes

