



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic test tube rack f. 2 x 12 tubes Ø 18 mm, 57 mm high

**€39.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Heidelberg Hospital model from our own workshop!

Sturdy, solid acrylic test tube racks with firmly glued side walls.

Heavy duty design with thick bottom.

Best material, original acrylic glass.

Hygienic, easy to wash.

Extremely durable, 1 year warranty on bonding and screwing.

Special sizes on request!

## Additional Information

No.	1-5104
Manufacturer (Brand)	neoLabLine
EAN	4058072035679
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	322 mm
Height	57 mm
Depth	64 mm
suitable for	Test tubes/tubes
Racks features	incubator suitable for different vessel types/volumes
Number of floors	1
DM Bore	18 mm
Squares (WxL)	2 x 12
Type frame	Test tube rack
for vessels with DM up to	18.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes
suitable for cooling in	Freezer/ freezer

