



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic test tube rack for 2 x 10 tubes 20 mm Ø, 57 mm high

**€25.02**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Heidelberg Hospital model from our own workshop!

Sturdy, solid acrylic test tube racks with firmly glued side walls.

Heavy duty design with thick bottom.

Best material, original acrylic glass.

Hygienic, easy to wash.

Extremely durable, 1 year warranty on bonding and screwing.

Special sizes on request!

## Additional Information

No.	1-5132
Manufacturer (Brand)	neoLabLine
EAN	4058072071462
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	312 mm
Height	57 mm
Depth	70 mm
suitable for	Test tubes/tubes
Racks features	incubator suitable for different vessel types/volumes
Number of floors	1
DM Bore	20 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 10
Type frame	Test tube rack
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

