



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




## Marienfeld test tube rack, 5 x 12 tubes with 16 mm Ø

**€9.00** As low as

**€7.65**

**plus VAT &  
Shipping**

### Product Images



## Description

---

**Test tube racks** made of white polypropylene with 3 shelves

- shatterproof
- suitable for water bath
- autoclavable
- stackable
- can be dismantled and stored in a space-saving way

## Additional Information

---

No.	MD-1421
Manufacturer (Brand)	Marienfeld
old neoLab article no.	E-0701
EAN	4058072156299
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	246 mm
Height	64 mm
Depth	104 mm
Racks features	stackable/collapsible suitable for incubator for different vessel types/volumes alphanumerically coded
DM Bore	16 mm
for vessels with volume up to	10 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	5 x 12
Type frame	Test tube rack
for vessels with DM up to	15.5 mm
for micro vessel/tube shape	for all shapes

