



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Marienfeld test tube rack, 4 x 10 tubes with 20 mm Ø

€9.00 As low as
€7.65
plus VAT &
Shipping

Product Images



Description

Test tube racks made of white polypropylene with 3 shelves

- shatterproof
- suitable for water bath
- autoclavable
- stackable
- can be dismantled and stored in a space-saving way

Additional Information

No.	MD-1422
Manufacturer (Brand)	Marienfeld
old neoLab article no.	E-0702
EAN	4058072160166
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	246 mm
Height	64 mm
Depth	104 mm
Racks features	stackable/collapsible suitable for incubator for different vessel types/volumes alphanumerically coded
DM Bore	20 mm
for vessels with volume up to	30 ml 50 ml
suitable for cooling in	Refrigerator freezer/ freezer
Squares (WxL)	4 x 10
Type frame	Test tube rack
for vessels with DM up to	19.5 mm
for micro vessel/tube shape	for all shapes

