



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® stainless steel test tube rack for 16 mm Ø, 4 x 12 rows

€61.10
plus VAT &
Shipping

Product Images



Description

Indestructible test tube racks for all purposes, made of 18/8 steel, polished.
For test tubes or centrifuge tubes.

Additional Information

No.	1-3131
Manufacturer (Brand)	neoLab
EAN	4058072035181
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Wide	250 mm
Height	100 mm
Depth	90 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
suitable for cooling in	Nitrogen tank
Tray size	17 x 17 mm
Squares (WxL)	4 x 12
Type frame	Test tube rack
for vessels with DM up to	17.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

