



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® stainless steel test tube rack for 12 mm Ø, 3 x 12 rows

**€49.50**  
plus VAT &  
Shipping

### Product Images



## Description

---

Indestructible test tube racks for all purposes, made of 18/8 steel, polished.  
For test tubes or centrifuge tubes.

## Additional Information

---

|                             |   |
|-----------------------------|---|
| No.                         | 1-3132  |
| Manufacturer (Brand)        | neoLab  |
| EAN                         | 4058072035198   |
| autoclavable at 121°C       | Yes   |
| Color                       | silver  |
| Material                    | Stainless steel   |
| Wide                        | 200 mm  |
| Height                      | 70 mm   |
| Depth                       | 60 mm   |
| suitable for                | Test tubes/tubes  |
| Racks features              | with special chemical resistance suitable for incubators for different vessel types/volumes |
| Number of floors            | 1   |
| suitable for cooling in     | Nitrogen tank   |
| Tray size                   | 14 x 14 mm  |
| Squares (WxL)               | 3 x 12  |
| Type frame                  | Test tube rack  |
| for vessels with DM up to   | 12.0 mm   |
| for micro vessel/tube shape | for all shapes  |
| for type (micro) vessel     | Test tubes/tubes  |

