



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® stainless steel test tube rack for 12 mm Ø, 4 x 12 rows

**€54.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Indestructible test tube racks for all purposes, made of 18/8 steel, polished.  
For test tubes or centrifuge tubes.

## Additional Information

---

No.	1-3133
Manufacturer (Brand)	neoLab
EAN	4058072035204
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Wide	200 mm
Height	70 mm
Depth	80 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
suitable for cooling in	Nitrogen tank
Tray size	14 x 14 mm
Squares (WxL)	4 x 12
Type frame	Test tube rack
for vessels with DM up to	12.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

