



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® test tube rack wood for 12 tubes on 2 shelves, hole Ø 20 mm

**€12.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

In line with the trend: ecological, renewable, produced in an environmentally friendly way!  
Test tube racks made of birch multiplex wood for 6 or 12 tubes.  
Waterproof glued (BFU100 according to DIN 68705 TI.3).

## Additional Information

---

No.	E-7053
Manufacturer (Brand)	neoLab
EAN	4058072014940
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	wood color
Material	Wood
Wide	192 mm
Height	184 mm
Depth	60 mm
suitable for	Test tubes/tubes
Racks features	with several tiers for different vessel types/volumes
Number of floors	2
DM Bore	20 mm
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Squares (WxL)	2 x 6
Type frame	Test tube rack
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

