



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® test tube rack wood for 6 tubes, hole Ø 20 mm, with 6 draining rods

€9.76
plus VAT &
Shipping

Product Images



Description

In line with the trend: ecological, renewable, environmentally friendly production!

Test tube racks made of birch multiplex wood for 6 or 12 glasses with 6 draining rods made of solid beech wood. Waterproof glued (BFU100 according to DIN 68705 Tl.3).

Additional Information

No.	E-7054
Manufacturer (Brand)	neoLab
EAN	4058072014957
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	wood color
Material	Wood
Wide	192 mm
Height	100 mm
Depth	65 mm
suitable for	Test tubes/tubes
Racks features	for different vessel types/volumes
Number of floors	1
DM Bore	20 mm
Squares (WxL)	1 x 6
Type frame	Test tube rack
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

