



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® test tube rack wood for 6 tubes, hole Ø 22 mm, with 6 draining rods

**€12.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

In line with the trend: ecological, renewable, environmentally friendly production!

Test tube racks made of birch multiplex wood for 6 or 12 glasses with 6 draining rods made of solid beech wood. Waterproof glued (BFU100 according to DIN 68705 TI.3).

## Additional Information

---

No.	E-7055
Manufacturer (Brand)	neoLab
EAN	4058072014964
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	wood color
Material	Wood
Wide	192 mm
Height	110 mm
Depth	80 mm
suitable for	Test tubes/tubes
Racks features	for different vessel types/volumes
Number of floors	1
DM Bore	22 mm
Squares (WxL)	1 x 6
Type frame	Test tube rack
for vessels with DM up to	22.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

