



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® test tube rack wood f. 12 tubes, 2 shelves, hole-Ø 20 mm, 6 draining rods

€13.70
plus VAT &
Shipping

Product Images



Description

In line with the trend: ecological, renewable, environmentally friendly production!

Test tube racks made of birch multiplex wood for 6 or 12 glasses with 6 draining rods made of solid beech wood. Waterproof glued (BFU100 according to DIN 68705 TI.3).

Additional Information

No.	E-7057
Manufacturer (Brand)	neoLab
EAN	4058072014988
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	wood color
Material	Wood
Wide	192 mm
Height	184 mm
Depth	80 mm
suitable for	Test tubes/tubes
Racks features	with several tiers for different vessel types/volumes
Number of floors	2
DM Bore	20 mm
Squares (WxL)	2 x 6
Type frame	Test tube rack
for vessels with DM up to	20.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

