



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® inclined rack for test tubes 13-16 mm Ø

**€77.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For easy and fast preparation of e.g. agar-media-slants.

On one side the tubes have an angle of 5° (useful for aerobic cultures).

On the other side the angle is 20° (well suited for non-aerobic cultures - reduces drying of substances during storage).

Substance-filled test tubes are autoclavable with the rack in the upright position.

40 spaces, alphanumerically labeled. Not stackable.

Dimensions: 267 x 152 x 130 mm (L x W x H).

Material: polypropylene,  
autoclavable up to 121°C

## Additional Information

No.	6-2320
Manufacturer (Brand)	neoLab
EAN	4058072025724
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	267 mm
Height	130 mm
Depth	152 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	16 mm
suitable for cooling in	Freezer/ freezer
Squares (WxL)	4 x 10
Type frame	Test tube rack
for vessels with DM up to	16.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Test tubes/tubes

