



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic Rack for 2 x 12 Eppendorf Reaction Tubes 2.2 ml

**€44.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Glued and screwed, Heidelberg Hospital model, heavy version. Height: 50 mm

## Additional Information

---

No.	1-5129
Manufacturer (Brand)	neoLabLine
EAN	4058072035808
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	50 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable incubator suitable
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	2 ml 2.2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	conical
for type (micro) vessel	Reaction vessels

