



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Acrylic Rack for 4 x 12 Eppendorf Reaction Tubes 2.2 ml

€39.36
plus VAT &
Shipping

Product Images



Description

Glued and screwed, Heidelberg Hospital model, heavy version. Height: 50 mm

Additional Information

No.	1-5130
Manufacturer (Brand)	neoLabLine
EAN	4058072035815
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	100 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable incubator suitable
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	2 ml 2.2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	4 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	conical
for type (micro) vessel	Reaction vessels

