



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic Racks 2 x12 hole, orange for Eppendorf Reaction Tubes

€46.90
plus VAT &
Shipping

Product Images



Description

Orange racks for 24 reaction tubes 1.5 ml.

Sturdy, solid acrylic racks with firmly glued and screwed side walls. Heavy duty design with thick base.

Side wall and top plate orange.

Best material, original acrylic glass

Hygienic, easy to wash.

Extremely durable.

Dimensions: 270 x 50 x 50 mm (L x W x H)

Additional Information

No.	1-5066
Manufacturer (Brand)	neoLabLine
EAN	4058072035617
autoclavable at 121°C	No
Color	orange
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	100 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable incubator suitable
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	1.5 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	4 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	conical
for type (micro) vessel	Reaction vessels

