



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Acrylic Racks blue for Eppendorf Reaction Tubes 2x6 hole

**€35.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Blue racks for 12 reaction tubes 1.5 ml.

Sturdy, solid acrylic racks with firmly glued and screwed side walls. Heavy duty design with thick bottom.

Best material, original acrylic glass.

Hygienic, easy to wash ,extremely durable.

Dimensions: 145 x 50 x 50 mm (L x W x H)

## Additional Information

---

|                               |   |
|-------------------------------|---|
| No.                           | 1-5061  |
| Manufacturer (Brand)          | neoLabLine                                    |
| EAN                           | 4058072035563                                 |
| autoclavable at 121°C         | No  |
| Color                         | blue  |
| Material                      | Polymethylmetacrylate (PMMA) (acrylic glass)  |
| Wide                          | 50 mm   |
| Height                        | 50 mm   |
| suitable for                  | Micro Vessels / Tubes                         |
| Racks features                | stackable/assemblable incubator suitable      |
| Number of floors              | 1   |
| DM Bore                       | 11 mm   |
| for vessels with volume up to | 1.5 ml  |
| suitable for cooling in       | Freezer/ freezer                              |
| Squares (WxL)                 | 2 x 6   |
| Type frame                    | Rack for reaction tubes. PCR tubes. Cryotubes |
| for vessels with DM up to     | 11.0 mm                                       |
| for micro vessel/tube shape   | conical                                       |
| for type (micro) vessel       | Reaction vessels                              |

