



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Acrylic Rack red for Eppendorf Reaction Tubes 2x6 hole

€35.10
plus VAT &
Shipping

Product Images



Description

Red racks for 12 reaction tubes 1.5 ml.

Stable, solid acrylic racks with firmly glued and screwed side walls. Heavy duty design with thick bottom.

Best material, original acrylic glass.

Hygienic, easy to wash ,extremely durable.

Dimensions: 145 x 50 x 50 mm (L x W x H)

Additional Information

No.	1-5062
Manufacturer (Brand)	neoLabLine
EAN	4058072035570
autoclavable at 121°C	No
Color	red
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	50 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable incubator suitable
Number of floors	1
DM Bore	11 mm
Capacity	12
Squares (WxL)	2 x 6
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for vessels with volume up to	1.5 ml
for micro vessel/tube shape	conical
for type (micro) vessel	Reaction vessels
suitable for cooling in	Freezer/ freezer

