



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Acrylic Racks 4 x12 hole, red for Eppendorf Reaction Tubes

**€95.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Red racks for 48 reaction tubes 1.5 ml.

Sturdy, solid acrylic racks with firmly glued and screwed side walls. Heavy duty design with thick bottom.

Side wall and top plate red.

Best material, original acrylic glass.

Hygienic, easy to wash.

Extremely durable.

Dimensions: 270 x 100 x 50 mm (L x W x H)

## Additional Information

No.	1-5065
Manufacturer (Brand)	neoLabLine
EAN	4058072035600
autoclavable at 121°C	No
Color	red
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	50 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable incubator suitable
Number of floors	1
DM Bore	11 mm
Squares (WxL)	2 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for vessels with volume up to	1.5 ml
for micro vessel/tube shape	conical
for type (micro) vessel	Reaction vessels
suitable for cooling in	Freezer/ freezer

