



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® acrylic rack for 5 ml Eppendorf tubes

€43.50
plus VAT &
Shipping

Product Images



Description

Sturdy, solid acrylic frame with firmly glued and screwed side parts.

- suitable for 5 ml Eppendorf reaction tubes
- 2 x 8 slots
- Holes are numerically marked
- easily washable
- extremely durable

Additional Information

No.	1-5125
Manufacturer (Brand)	neoLabLine
EAN	4058072035792
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Wide	50 mm
Height	55 mm
Number of floors	1
DM Bore	17 mm
for vessels with volume up to	5 ml
suitable for cooling in	Refrigerator
Squares (WxL)	2 x 8
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	16.0 mm
for micro vessel/tube shape	conical

