



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand reaction tube stand, gray for 20 reaction tubes, 1.5 ml

€11.95
plus VAT &
Shipping

Product Images



Description

Gray reaction tube stand made of PP. Positions (numbered) for 20 reaction tubes with 1.5 ml. Autoclavable at 121°C.

Additional Information

No.	BN-2151
Manufacturer (Brand)	Brand
old neoLab article no.	1-6819
EAN	4058072213756
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	70 mm
Height	37 mm
Racks features	alphanumeric coded
DM Bore	11.5 mm
for vessels with volume up to	1.5 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Squares (WxL)	2 x 10
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	conical

