



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Beta safety rack for 4 x 6 reaction tubes 1.5/2.0 ml

€197.40
plus VAT &
Shipping

Product Images



Description

Solid cast acrylic block with holes to hold 4 x 6 reaction tubes each.
Beta-radiation absorbing with hold-up lid and **warning label**

Order No. 4-2006

For reaction vessels 1.5/2.0 ml

Additional Information

No.	4-2006
Manufacturer (Brand)	neoLab
EAN	4058072022723
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Wide	110 mm
Height	150 mm
Racks features	with closure/lid with beta rays protection
Number of floors	1
DM Bore	11.5 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Refrigerator
Squares (WxL)	4 x 6
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	conical

