



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Combi-Rack for PCR tubes and reaction vessels with lid, yellow


Special Price

€26.50 was

€33.40

**plus VAT &
Shipping**

Product Images



Description

Stable rack made of polypropylene for storage of different sized tubes, with transparent PP lid.

24 places for reaction tubes 1.5 ml ,16 places for reaction tubes 0.5 ml.

Removable insert for 96 (8 x 12) PCR tubes.

Stackable, autoclavable up to 121°C.

Dimensions:

210 x 115 mm (L x W)

120 x 85 mm (L x W) PCR insert

Additional Information

No.	2-2760
Manufacturer (Brand)	neoLab
EAN	4058072068912
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
sterile	No
Wide	115 mm
Height	78 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance with closure/lid suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
Squares (WxL)	variable
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels
suitable for cooling in	Ultra Deep Freezer / Deep Freezer

