



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Combi-Rack for PCR tubes and reaction vessels, with lid, pink


Special Price

**€26.50** was

**€33.40**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Stable rack made of polypropylene for storage of different sized tubes, with transparent PP lid.

24 places for reaction tubes 1.5 ml ,16 places for reaction tubes 0.5 ml.

Removable insert for 96 (8 x 12) PCR tubes.

Stackable, autoclavable up to 121°C.

Dimensions:

210 x 115 mm (L x W)

120 x 85 mm (L x W) PCR insert

## Additional Information

No.	2-2763
Manufacturer (Brand)	neoLab
EAN	4058072068943
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	pink
Material	Polypropylene (PP)
Wide	115 mm
Height	78 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance with closure/lid suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
Squares (WxL)	variable
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels
suitable for cooling in	Ultra Deep Freezer / Deep Freezer

