



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Combi Rack for reaction tubes, 96 places, yellow

**€29.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Double-sided rack for reaction tubes**

**Top side** for 1.5 and 2 ml tubes

**Bottom side** for 0.5 ml tubes

**Dimensions:** 210 x 115 x 28 mm

**Material** : PP, autoclavable up to 121°C

Accessories: Transparent lid made of PP (Order No. 2-2595)

## Additional Information

No.	2-2570
Manufacturer (Brand)	neoLabLine
EAN	4058072061333
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	115 mm
Height	28 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/collapsible with special chemical resistance suitable for incubator for different vessel types/volumes closure/lid available as accessory
Number of floors	1
Squares (WxL)	8 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels
suitable for cooling in	Freezer/ freezer

