



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Combi Rack for reaction tubes, 96 places, transparent

€29.50
plus VAT &
Shipping

Product Images



Description

Double-sided rack for reaction tubes

Top side for 1.5 and 2 ml tubes

Bottom side for 0.5 ml tubes

Dimensions: 210 x 115 x 28 mm

Material : PP, autoclavable up to 121°C

Accessories: Transparent lid made of PP (Order No. 2-2595)

Additional Information

No.	2-2576
Manufacturer (Brand)	neoLabLine
EAN	4058072068165
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Wide	115 mm
Height	28 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/collapsible with special chemical resistance suitable for incubator for different vessel types/volumes closure/lid available as accessory
Number of floors	1
Capacity	96
Squares (WxL)	8 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels
suitable for cooling in	Freezer/ freezer

