



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Combi Rack with snap-on lid for 96 reaction tubes, yellow

**€15.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Storage rack with snap-on lid for

- 96 reaction tubes 0.5 ml on one side or
- 96 reaction tubes 1.5/2.0 ml on the other side.

A grid on the side for the 0.5 ml tubes makes it easier to keep track.

The side for 1.5/2.0 ml tubes has raised hole borders to make it easier to lift the tubes out of the rack.

Material: Polypropylene

Spaces: 8 x 12

Dimensions: 246 x 121 x 50 mm

## Additional Information

No.	2-1837
Manufacturer (Brand)	neoLab
EAN	4058072058173
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	121 mm
Height	50 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance with closure/lid suitable for incubator for different vessel types/volumes
Number of floors	1
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	8 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels

