



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Combi Rack for PCR Tubes, PP, 24 x 0.2 ml, 14 x 0.5 ml, 12 x 1.5/2.0 ml, pink

Special Price

€5.00 was

€7.50

plus VAT &
Shipping

Product Images



Description

Combi rack for PCR tubes

The pink three-tier combirack for PCR tubes is made of polypropylene and therefore autoclavable. Due to the practical three-tiered shape, the Kombirack offers space for the following reaction tubes:

- 24 x 0.2 ml
- 14 x 0.5 ml
- 12 x 1.5/2.0 ml

Additional Information

No.	7-3097
Manufacturer (Brand)	neoLab
EAN	4058072087463
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	pink
Material	Polypropylene (PP)
Racks features	stackable/connectable with several tiers for different vessel types/volumes
Number of floors	3
Capacity	50
Squares (WxL)	2 x 12 2 x 7 2 x 6
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	13.0 mm
for vessels with volume up to	0.2 ml 0.5 ml 1.5 ml 2 ml
for micro vessel/tube shape	conical round
suitable for cooling in	Refrigerator

