



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene rack made of PC for cryo tubes, 5 x 10 places

€49.20 As low
as
€39.75
plus VAT &
Shipping

Product Images



Description

Cryo-tube rack

For 5 x 10 cryo-tubes 1.2 and 2 ml, with marking for identification of the respective tube, with peg in the opening that holds the cryo-tube in place.

Autoclavable, made of polycarbonate.

Additional Information

No.	CD-1536
Manufacturer (Brand)	Nalgene
old neoLab article no.	3-6002
EAN	4058072219406
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polycarbonate (PC)
Wide	102 mm
Height	197 mm
suitable for	Micro Vessels / Tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	13.5 mm
Squares (WxL)	5 x 10
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	13.0 mm
for vessels with volume up to	1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Freezer/ freezer

