



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo tube holder aluminum, 290 mm long, 12 pcs./pack

€32.70
plus VAT &
Shipping

Product Images



Description

Aluminum holder for clamping the cryovials, up to 5 pieces of 1.2 ml or 2 ml cryotubes with 12.5 mm Ø.
For storing the tubes in liquid nitrogen tanks or Dewar flasks.
Total length: 290 mm

Additional Information

No.	7-8029
Manufacturer (Brand)	neolab
EAN	4058072041847
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Length	290 mm
Racks features	for different vessel types/volumes
for vessels with volume up to	2 ml
suitable for cooling in	Nitrogen tank
Squares (WxL)	5 x 1
Type frame	Cryotube rack
for vessels with DM up to	12.5 mm

