



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® cryo tube workstation 5 x 10 places,  
red, 4 pcs./pack**

**€83.90  
plus VAT &  
Shipping**

## Product Images



## Description

---

Handy, autoclavable rack for cryotubes.

Small pegs in the recesses hold the tubes in place when opening or closing them.

Small rubber feet prevent slipping.

Molded handles on both sides facilitate safe transport.

Material: PP autoclavable

Openings: 5 rows with 10 places each

Dimensions: 20 x 10 x 2.5 cm high

## Additional Information

---

No.	7-8025
Manufacturer (Brand)	neoLab
EAN	4058072041823
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Wide	100 mm
Height	25 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	12.5 mm
for vessels with volume up to	5 ml
suitable for cooling in	Refrigerator
Squares (WxL)	5 x 10
Type frame	Cryotube rack
for vessels with DM up to	12.5 mm
for micro vessel/tube shape	for all shapes

