



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cryo tube workstation 5 x 10 places,
pink, 4 pcs./pack**

€83.90
plus VAT &
Shipping

Product Images

Description

Handy, autoclavable rack for cryotubes.

Small pegs in the recesses hold the tubes in place when opening or closing them.

Small rubber feet prevent slipping.

Molded handles on both sides facilitate safe transport.

Material: PP autoclavable

Openings: 5 rows with 10 places each

Dimensions: 20 x 10 x 2.5 cm high

Additional Information

No.	7-8027
Manufacturer (Brand)	neoLab
EAN	4058072041830
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	pink
Material	Polypropylene (PP)
Wide	100 mm
Height	25 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	12.5 mm
Squares (WxL)	5 x 10
Type frame	Cryotube rack
for vessels with DM up to	12.5 mm
for vessels with volume up to	5 ml
for micro vessel/tube shape	for all shapes
suitable for cooling in	Refrigerator

