



Thermo Scientific Matrix™ Reagent Reservoirs

Product Images



Description

Matrix™ Reagent Reservoirs Create space for a variety of applications with the proprietary, innovative reagent reservoirs. Packaging: 10 bags each 10 capacity (metric): 100 ml Type: Reagent Reservoir Unit Size: Pack of 100

Additional Information

No.	TA-0404
Manufacturer (Brand)	Thermo Scientific
EAN	4066292338852
Gross weight	1.484
Color	Weiß
Material	Polystyrol (PS)
Country of origin	USA
sterile	No
Packaging volume	0.01300000 ccm
Customs tariff number	39269097

