



Thermo Scientific Matrix™ Reagent Reservoirs

Product Images



Description

Matrix™ Reagent Reservoirs, 25 ml, with dividers, pouches Create space for a variety of applications with the specially developed, innovative reagent reservoirs. Packaging: 10 bags of 10 capacity each (metric): 25 ml Type: Subdivided Reagent Reservoirs Unit Size: Pack of 100

Additional Information

No.	TA-0408
Manufacturer (Brand)	Thermo Scientific
EAN	4066292338883
Gross weight	1
Color	Weiß
Material	Polystyrol (PS)
Country of origin	USA
sterile	No
Packaging volume	0.00900000 ccm
Customs tariff number	39269097

