



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® Rack for reaction tubes (1.5/2.0 ml), 24 places, yellow, PC, 1 pc

€4.16
plus VAT &
Shipping

Product Images



Description

This yellow rack for 1.5 and 2.0 ml reaction tubes features an ergonomic design. With the raised bottom, the reaction tubes can be easily removed. The rack is alphanumerically coded and can be stacked (empty or filled).

The rack is made of polycarbonate and is autoclavable.

Additional Information

No.	4-0019
Manufacturer (Brand)	Moonlab
EAN	4058072085162
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polycarbonate (PC)
Wide	100 mm
Height	59 mm
Racks features	for different vessel types/volumes
Number of floors	1
DM Bore	13 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Refrigerator
Squares (WxL)	4 x 6
for vessels with DM up to	13.0 mm
for micro vessel/tube shape	conical round

