



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® Reversible rack for reaction tubes with lids (1.5/2.0/0.5 ml), 96 places (12x8), yellow, PP, 1 pc, 219x150x50 mm**

**€5.49  
plus VAT &  
Shipping**

## Product Images



## Description

---

Simple and stable, double-sided racks with transparent lid for 96 +6 reaction tubes.

**Top side:** 12 x 8 rows for 1.5 or 2.0 ml reaction vessels ( 96 places) with additional 6 places for 0.5 ml reaction vessels.

**Bottom side:** 12 x 8 rows for 0.5 ml reaction vessels ( 96 places) with additional 6 places for 1.5/2.0 ml

## Additional Information

---

No.	4-0010
Manufacturer (Brand)	Moonlab
EAN	4058072085711
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
Material	Polypropylene (PP)
Wide	150 mm
Height	50 mm
Racks features	with closure/lid for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	13.5 mm
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Refrigerator
Squares (WxL)	8 x 12
for vessels with DM up to	13.0 mm 9.0 mm
for micro vessel/tube shape	conical round

