



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® reaction tube rack 4 x 6 places, 4 pcs./pack

€50.60
plus VAT &
Shipping

Product Images



Description

Inexpensive stand for reaction tubes 1.5 ml capacity or equivalent.

Floatable.

Made of polypropylene, autoclavable up to 121°C.

With handle for removal from water bath. Suspension rack.

4 x 6 places (24) with 10 mm Ø

Dimensions: 120 x 80 mm

Additional Information

No.	1-3114
Manufacturer (Brand)	neolab
EAN	4058072035075
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Wide	80 mm
suitable for	Micro Vessels / Tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	4 x 6
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels

