




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rota-Rack® Duo for PCR and reaction tubes

Special Price

€48.70 was

€61.30

plus VAT &
Shipping

Product Images



Description

Small in format, big in combining!

Ideal for working with different vessel sizes side by side, e.g. when diluting or concentrating.

Two rotating, interconnected racks made of PP provide space for 96 vessels with 0.2 ml and 48 places each for vessels with 0.5 and 1.5 ml contents.

The 0.2 ml wells have the standard format for 96-well plates.

Quick identification by alphanumeric marking.

Autoclavable, temperature resistant from 0 to +121°C.

Color: blue and purple

Dimensions: 298 x 118 x 34 mm (L x D x H)

Additional Information

No.	2-6626
Manufacturer (Brand)	neoLab
EAN	4058072008222
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue violet
Material	Polypropylene (PP)
Wide	118 mm
Height	34 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
for vessels with volume up to	0.2 ml 0.5 ml 1.5 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels

