



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand Rack with grid for tubes 1.2 ml, empty, non-sterile

**€68.75**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to  $-80^{\circ}\text{C}$ .

The lids made of PE are not suitable for autoclaving.

Dimensions of the tubes:

- Outer diameter: 8.8 mm

- Height: 45 mm

## Additional Information

---

|                                       |                    |
|---------------------------------------|--------------------|
| No.                                   | BN-2287            |
| Manufacturer (Brand)                  | Brand              |
| old neoLab article no.                | 7-5244             |
| EAN                                   | 4058072214968      |
| Transport temperature                 | Room temperature   |
| autoclavable at $121^{\circ}\text{C}$ | Yes                |
| Color                                 | White              |
| Material                              | Polypropylene (PP) |
| Height                                | 50 mm              |
| suitable for                          | Reaction vessels   |
| Racks features                        | with closure/lid   |

