



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® slot rack, 9 slots à 25.4 mm width

€87.70
plus VAT &
Shipping

Product Images



Description

Slotted rack

For storing or cold storage of small vials of blood serum, chemical reagents, enzyme reagents, etc. in the refrigerator. Can also be used for sorting pipettes, stirrers or similar in drawers.

Vacuum formed.

Dimensions: 330 x 330 x 44 mm

Additional Information

No.	6-2240
Manufacturer (Brand)	neoLab
EAN	4058072025564
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Wide	330 mm
Height	44 mm
Length	330 mm
Depth	330 mm
suitable for	Vials/ Vials
Number of compartments	9
for vessels with DM up to	25.0 mm
for type vial/vial	Vial. 25 mm Ø
suitable for cooling in	Freezer/ freezer

