



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® slot rack, 8 rows of 32 mm width each

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

### Slotted rack

For storing or cold storage of small vials of blood serum, chemical reagents, enzyme reagents, etc. in the refrigerator. Can also be used for sorting pipettes, stirrers or similar in drawers.

Vacuum formed.

Dimensions: 330 x 330 x 44 mm

## Additional Information

---

No.	6-2241
Manufacturer (Brand)	neoLab
EAN	4058072025571
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Wide	330 mm
Height	44 mm
Length	330 mm
Depth	330 mm
suitable for	Vials/ Vials
Number of compartments	8
suitable for cooling in	Freezer/ freezer
for vessels with DM up to	30.0 mm
for type vial/vial	Vial. 30 mm Ø

