



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Reaction Tube Rack, yellow, PP, 2 x 10 tubes 1.5/2.0 ml

**€9.76** As low as  
**€8.22**  
plus VAT &  
Shipping

### Product Images



## Description

---

The neoLab reaction tube stand is practical and chic at the same time. Available in different colors (blue, yellow, green, red) and made of durable polypropylene, the neoLab reaction tube stand is temperature resistant from  $-90^{\circ}\text{C}$  to  $+120^{\circ}\text{C}$  and autoclavable. The feet of the reaction vessel rack are secured with knobs against slipping on the tabletop, but should be removed for autoclaving. With a division of 2 x 10 rows, the neoLab reaction tube rack offers space for 20 x 1.5 or 2.0 ml reaction tubes or eppis. Conical recesses guarantee a firm hold of the reaction vessels in the rack. The wedge shape also makes it easy to reach and process the reaction tubes in the back row. All openings are numbered.

## Additional Information

No.	2-1969
Manufacturer (Brand)	neoLabLine
EAN	4058072059132
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	50 mm
Height	73 mm
suitable for	Micro Vessels / Tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 10
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels

