



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Top Rack for 80 reaction tubes 1.5 ml, yellow

€6.32
plus VAT &
Shipping

Product Images



Description

Yellow storage rack with 80 numerically labeled spaces for reaction tubes 1.5 and 2.0 ml.

With handle strips; in various colors.

Temperature resistant from -80°C to +125°C.

Material: polypropylene, autoclavable

Additional Information

| | |
|-------------------------------|--|
| No. | 2-2560 |
| Manufacturer (Brand) | neolabLine |
| EAN | 4058072061258 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | yellow |
| Material | Polypropylene (PP) |
| Wide | 65 mm |
| Height | 28 mm |
| suitable for | Micro Vessels / Tubes |
| Racks features | stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes |
| Number of floors | 1 |
| DM Bore | 11 mm |
| for vessels with volume up to | 1.5 ml 2 ml |
| suitable for cooling in | Ultra Deep Freezer / Deep Freezer |
| Squares (WxL) | 5 x 16 |
| Type frame | Rack for reaction tubes. PCR tubes. Cryotubes |
| for vessels with DM up to | 11.0 mm |
| for micro vessel/tube shape | for all shapes |
| for type (micro) vessel | Reaction vessels |

