



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Top Rack for 80 reaction tubes 1.5 ml, violet

**€6.32**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Purple storage rack with 80 numerically labeled spaces for 1.5 and 2.0 ml reaction tubes.

With handle strips; in various colors.

Temperature resistant from -80°C to +125°C.

Material: polypropylene, autoclavable

## Additional Information

---

No.	2-2563
Manufacturer (Brand)	neolabLine
EAN	4058072061289
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	purple
Material	Polypropylene (PP)
Wide	65 mm
Height	28 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	11 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Ultra Deep Freezer / Deep Freezer
Squares (WxL)	5 x 16
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels

