



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Top Rack for 80 reaction tubes 1.5 ml, transparent

€6.32
plus VAT &
Shipping

Product Images



Description

Transparent storage rack with 80 numerically labeled spaces for 1.5 and 2.0 ml reaction tubes.

With handle strips; in various colors.

Temperature resistant from -80°C to +125°C.

Material: polypropylene, autoclavable

Additional Information

No.	2-2566
Manufacturer (Brand)	neoLabLine
EAN	4058072061319
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Wide	65 mm
Height	28 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	11 mm
Squares (WxL)	5 x 16
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	11.0 mm
for vessels with volume up to	1.5 ml 2 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels
suitable for cooling in	Ultra Deep Freezer / Deep Freezer

