




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Work-2-Store™ Storage Box for 1.5/20, ml

€38.00
plus VAT &
Shipping

Product Images



Description

With this cryobox you can store and process samples at the same time. The patented opening and closing mechanism makes it easy to expand the working area. After sample processing, simply push the rows back together and place the cryobox back in the freezer.

- immediately more space to work
- Avoidance of cross contamination
- opened reaction vessels do not obstruct each other
- stackable
- Alphanumeric coding
- Temperature range from -80 to 121°C (autoclavable)
- Dimensions closed with lid (LxWxH): 133x133x53 mm
- Dimensions opened (LxWxH): 310x128x335 mm

Capacities:

Order No. 2-6632

For 1.5 ml tubes: 64 places

For 2.0 ml tubes: 64 places

Additional Information

No.	2-6632
Manufacturer (Brand)	neoLab
EAN	4058072008284
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Wide	133 mm
Height	53 mm
Depth	133 mm
Racks features	stackable/assemblable with closure/lid alphanumerically coded
Number of floors	1
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator
Squares (WxL)	8 x 8
Type frame	Cryotube rack
for vessels with DM up to	11.0 mm

